

LensHOOKE®

R10

Sperm DNA Fragmentation
Test Kit (SCD Assay)



Certificate
FDA



Sperm DNA Fragmentation




Save Time

LensHooke® R10 develops the new in-gel denature technology to reduce 40% test time. It could allow you to complete one sample within 40 minutes.


In-Gel Denature Technology

Patented
TWI775252



Clear Image

LensHooke® R10 provides you with a clear image to reduce the subjective issue by manual.



Content Includes (10 Test / Box)

- **DNA Denaturant Reagent**
Hydrochloric acid, 2 mL
- **Lysis Solution**
Reducing agent, 15 mL
- **Stain Solution A**
Wright-Giemsa dye, 14 mL
- **Stain Solution B**
Phosphate, 14 mL
- **Microcentrifuge Tubes with Agarose Gel**
Low melting point agarose gel, 50 µL
- **Microcentrifuge Tubes with Agarose Gel (for Quality Control)**
Low melting point agarose gel, 0.3 mL
- **Pretreated glass slide**
Normal-melting agarose gel
- **Float**
Φ 6 cm

Product Specification

- **Operation Time**
40 minutes
- **Operation Environment**
15°C~38°C (59~100°F)
- **Storage Environment**
2°C~8°C (36°F~46°F)
- **Shelf Life**
Sealed for 24 months
Opened for 9 months

